



TESPA1 blocking peptide (CDBP6264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TESPA1 thymocyte expressed, positive selection associated 1 [Homo sapiens (human)]
Official Symbol	TESPA1
Synonyms	TESPA1; thymocyte expressed, positive selection associated 1; HSPC257; KIAA0748; protein TESPA1; thymocyte-expressed positive selection-associated protein 1
Entrez Gene ID	9840
mRNA Refseq	NM_001098815
Protein Refseq	NP_001092285
UniProt ID	A2RU30
Function	receptor binding