



THAP11 blocking peptide (CDBP6271)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene contains a THAP domain, which is a conserved DNA-binding domain that has striking similarity to the site-specific DNA-binding domain (DBD) of Drosophila P element transposases. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	THAP11 THAP domain containing 11 [Homo sapiens (human)]
Official Symbol	THAP11
Synonyms	THAP11; THAP domain containing 11; RONIN; CTG-B43a; CTG-B45d; HRIHFB2206; THAP domain-containing protein 11
Entrez Gene ID	57215
mRNA Refseq	NM_020457
Protein Refseq	NP_065190

UniProt ID	Q96EK4
Function	DNA binding; protein binding; zinc ion binding