



TIAF1 blocking peptide (CDBP6275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TIAF1 TGFB1-induced anti-apoptotic factor 1 [Homo sapiens (human)]
Official Symbol	TIAF1
Synonyms	TIAF1; TGFB1-induced anti-apoptotic factor 1; MAJN; SPR210; TGF-beta-1-induced antiapoptotic factor 1; molecule associated with Jak-3 N-terminal; 12 kDa TGF-beta-1-induced antiapoptotic factor
Entrez Gene ID	9220
mRNA Refseq	NM_004740
Protein Refseq	NP_004731
UniProt ID	O95411

Functionmolecular_function
