



TIMD4 blocking peptide (CDBP6279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TIMD4 T-cell immunoglobulin and mucin domain containing 4 [Homo sapiens (human)]
Official Symbol	TIMD4
Synonyms	TIMD4; T-cell immunoglobulin and mucin domain containing 4; TIM4; SMUCKLER; T-cell immunoglobulin and mucin domain-containing protein 4; TIM-4; TIMD-4; T-cell membrane protein 4; T-cell immunoglobulin mucin receptor 4
Entrez Gene ID	91937
mRNA Refseq	NM_001146726
Protein Refseq	NP_001140198
UniProt ID	Q96H15