



TICAM2 blocking peptide (CDBP6284)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TIRP is a Toll/interleukin-1 receptor (IL1R; MIM 147810) (TIR) domain-containing adaptor protein involved in Toll receptor signaling (see TLR4; MIM 603030).[supplied by OMIM, Apr 2004]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TICAM2 toll-like receptor adaptor molecule 2 [Homo sapiens (human)]
Official Symbol	TICAM2
Synonyms	TICAM2; toll-like receptor adaptor molecule 2; TIRP; TRAM; TIRAP3; MyD88-4; TICAM-2; TIR domain-containing adapter molecule 2; cytoplasmic adaptor; TRIF-related adaptor molecule; NF-kappa-B-activating protein 502; toll-like receptor adaptor protein 3; putative NF-kappa-B-activating protein 502; toll/interleukin-1 receptor (TIR) domain-containing adapter protein
Entrez Gene ID	353376
mRNA Refseq	NM_021649

Protein Refseq	NP_067681
UniProt ID	Q86XR7
Pathway	Activated TLR4 signalling; Activation of IRF3/IRF7 mediated by TBK1/IKK epsilon; Hepatitis B; IKK complex recruitment mediated by RIP1; IL1-mediated signaling events; Immune System; Innate Immune System; MyD88-independent cascade
Function	signal transducer activity