



Tlr11 blocking peptide (CDBP6289)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Tlr11 toll-like receptor 11 [Mus musculus (house mouse)]
Official Symbol	Tlr11
Synonyms	Tlr11; toll-like receptor 11; Gm287
Entrez Gene ID	239081
mRNA Refseq	NM_205819
Protein Refseq	NP_991388
UniProt ID	Q6R5P0
Pathway	Toxoplasmosis

Function

protein binding; transmembrane signaling receptor activity
