



TMEM106A blocking peptide (CDBP6307)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM106A transmembrane protein 106A [Homo sapiens (human)]
Official Symbol	TMEM106A
Synonyms	TMEM106A; transmembrane protein 106A
Entrez Gene ID	113277
mRNA Refseq	NM_001291586
Protein Refseq	NP_001278515
UniProt ID	Q96A25
Function	molecular_function