



TMEM106B blocking peptide (CDBP6308)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM106B transmembrane protein 106B [Homo sapiens (human)]
Official Symbol	TMEM106B
Synonyms	TMEM106B; transmembrane protein 106B
Entrez Gene ID	54664
mRNA Refseq	NM_001134232
Protein Refseq	NP_001127704
UniProt ID	Q9NUM4
Function	molecular_function; protein binding