



TMEM107 blocking peptide (CDBP6309)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM107 transmembrane protein 107 [Homo sapiens (human)]
Official Symbol	TMEM107
Synonyms	TMEM107; transmembrane protein 107; GRVS638; PRO1268
Entrez Gene ID	84314
mRNA Refseq	NM_032354
Protein Refseq	NP_115730
UniProt ID	Q6UX40
Function	molecular_function