



TMEM135 blocking peptide (CDBP6310)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM135 transmembrane protein 135 [Homo sapiens (human)]
Official Symbol	TMEM135
Synonyms	TMEM135; transmembrane protein 135; PMP52; peroxisomal membrane protein 52
Entrez Gene ID	65084
mRNA Refseq	NM_001168724
Protein Refseq	NP_001162195
UniProt ID	Q86UB9