



TMEM147 blocking peptide (CDBP6311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM147 transmembrane protein 147 [Homo sapiens (human)]
Official Symbol	TMEM147
Synonyms	TMEM147; transmembrane protein 147; NIFIE14; protein NIFIE 14; seven transmembrane domain protein
Entrez Gene ID	10430
mRNA Refseq	NM_001242597
Protein Refseq	NP_001229526
UniProt ID	Q9BVK8