



TMEM18 blocking peptide (CDBP6314)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM18 transmembrane protein 18 [Homo sapiens (human)]
Official Symbol	TMEM18
Synonyms	TMEM18; transmembrane protein 18
Entrez Gene ID	129787
mRNA Refseq	NM_152834
Protein Refseq	NP_690047
UniProt ID	Q96B42
Function	DNA binding