



TMEM184A blocking peptide (CDBP6315)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM184A transmembrane protein 184A [Homo sapiens (human)]
Official Symbol	TMEM184A
Synonyms	TMEM184A; transmembrane protein 184A
Entrez Gene ID	202915
mRNA Refseq	NM_001097620
Protein Refseq	NP_001091089
UniProt ID	Q6ZMB5