



TMEM184B blocking peptide (CDBP6316)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM184B transmembrane protein 184B [Homo sapiens (human)]
Official Symbol	TMEM184B
Synonyms	TMEM184B; transmembrane protein 184B; FM08; HS5O6A; C22orf5; putative MAPK-activating protein FM08
Entrez Gene ID	25829
mRNA Refseq	NM_001195071
Protein Refseq	NP_001182000
UniProt ID	Q9Y519