



TMEM184C blocking peptide (CDBP6317)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM184C transmembrane protein 184C [Homo sapiens (human)]
Official Symbol	TMEM184C
Synonyms	TMEM184C; transmembrane protein 184C; TMEM34; transmembrane protein 34
Entrez Gene ID	55751
mRNA Refseq	NM_018241
Protein Refseq	NP_060711
UniProt ID	Q9NVA4