



TMEM192 blocking peptide (CDBP6318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM192 transmembrane protein 192 [Homo sapiens (human)]
Official Symbol	TMEM192
Synonyms	TMEM192; transmembrane protein 192
Entrez Gene ID	201931
mRNA Refseq	NM_001100389
Protein Refseq	NP_001093859
UniProt ID	Q8IY95
Function	protein homodimerization activity