



TMEM204 blocking peptide (CDBP6319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C16ORF30 plays a role in cell adhesion and cellular permeability at adherens junctions (Kearsey et al., 2004 [PubMed 15206924]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM204 transmembrane protein 204 [Homo sapiens (human)]
Official Symbol	TMEM204
Synonyms	TMEM204; transmembrane protein 204; CLP24; C16orf30; claudin-like protein 24
Entrez Gene ID	79652
mRNA Refseq	NM_001256541
Protein Refseq	NP_001243470
UniProt ID	Q9BSN7