



TMEM214 blocking peptide (CDBP6321)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM214 transmembrane protein 214 [Homo sapiens (human)]
Official Symbol	TMEM214
Synonyms	TMEM214; transmembrane protein 214
Entrez Gene ID	54867
mRNA Refseq	NM_001083590
Protein Refseq	NP_001077059
UniProt ID	Q6NUQ4