



TMEM38A blocking peptide (CDBP6322)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM38A transmembrane protein 38A [Homo sapiens (human)]
Official Symbol	TMEM38A
Synonyms	TMEM38A; transmembrane protein 38A; TRICA; TRIC-A; trimeric intracellular cation channel type A
Entrez Gene ID	79041
mRNA Refseq	NM_024074
Protein Refseq	NP_076979
UniProt ID	Q9H6F2
Function	potassium channel activity