



TMEM38B blocking peptide (CDBP6323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an intracellular monovalent cation channel that functions in maintenance of intracellular calcium release. Mutations in this gene may be associated with autosomal recessive osteogenesis. [provided by RefSeq, Oct 2012]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM38B transmembrane protein 38B [Homo sapiens (human)]
Official Symbol	TMEM38B
Synonyms	TMEM38B; transmembrane protein 38B; OI14; TRICB; TRIC-B; C9orf87; D4Ert89e; bA219P18.1; trimeric intracellular cation channel type B
Entrez Gene ID	55151
mRNA Refseq	NM_018112
Protein Refseq	NP_060582

UniProt ID	Q9NVV0
Function	potassium channel activity