



TMEM59 blocking peptide (CDBP6324)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM59 transmembrane protein 59 [Homo sapiens (human)]
Official Symbol	TMEM59
Synonyms	TMEM59; transmembrane protein 59; C1orf8; HSPC001; liver membrane-bound protein
Entrez Gene ID	9528
mRNA Refseq	NM_004872
Protein Refseq	NP_004863
UniProt ID	Q9BXS4
Function	endopeptidase activity; protein binding