



TMEM59L blocking peptide (CDBP6325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a predicted type-I membrane glycoprotein. The encoded protein may play a role in functioning of the central nervous system. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMEM59L transmembrane protein 59-like [Homo sapiens (human)]
Official Symbol	TMEM59L
Synonyms	TMEM59L; transmembrane protein 59-like; BSMAP; C19orf4; brain-specific membrane-anchored protein
Entrez Gene ID	25789
mRNA Refseq	NM_012109
Protein Refseq	NP_036241

