



# TMEM88 blocking peptide (CDBP6327)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                          |
|-----------------|--------------------------------------------------------------------------|
| Gene Name       | <a href="#">TMEM88 transmembrane protein 88 [ Homo sapiens (human) ]</a> |
| Official Symbol | TMEM88                                                                   |
| Synonyms        | TMEM88; transmembrane protein 88                                         |
| Entrez Gene ID  | <a href="#">92162</a>                                                    |
| mRNA Refseq     | <a href="#">NM_203411</a>                                                |
| Protein Refseq  | <a href="#">NP_981956</a>                                                |
| UniProt ID      | Q6PEY1                                                                   |