



TMIGD1 blocking peptide (CDBP6328)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMIGD1 transmembrane and immunoglobulin domain containing 1 [Homo sapiens (human)]
Official Symbol	TMIGD1
Synonyms	TMIGD1; transmembrane and immunoglobulin domain containing 1; TMIGD; UNQ9372; transmembrane and immunoglobulin domain-containing protein 1; AWKS9372
Entrez Gene ID	388364
mRNA Refseq	NM_206832
Protein Refseq	NP_996663
UniProt ID	Q6UXZ0