



TMIGD2 blocking peptide (CDBP6329)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMIGD2 transmembrane and immunoglobulin domain containing 2 [Homo sapiens (human)]
Official Symbol	TMIGD2
Synonyms	TMIGD2; transmembrane and immunoglobulin domain containing 2; IGPR1; IGPR-1; transmembrane and immunoglobulin domain-containing protein 2; CD28 homolog; immunoglobulin and proline-rich receptor 1; immunoglobulin-containing and proline-rich receptor 1; transmembrane and immunoglobulin domain-containing protein 2 variant 2; transmembrane and immunoglobulin domain-containing protein 2 variant 3
Entrez Gene ID	126259
mRNA Refseq	NM_001169126
Protein Refseq	NP_001162597

UniProt ID	Q96BF3
Function	coreceptor activity; protein binding