



TMPRSS11A blocking peptide (CDBP6332)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TMPRSS11A transmembrane protease, serine 11A [Homo sapiens (human)]
Official Symbol	TMPRSS11A
Synonyms	TMPRSS11A; transmembrane protease, serine 11A; ECRG1; transmembrane protease serine 11A; airway trypsin-like protease 1; epidermal type II transmembrane serine protease; epidermal type-II transmembrane serine protease; esophageal cancer-susceptibility gene 1 protein
Entrez Gene ID	339967
mRNA Refseq	NM_001114387
Protein Refseq	NP_001107859
UniProt ID	Q6ZMR5

Functionserine-type endopeptidase activity
