



TOMM70A blocking peptide (CDBP6344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes an import receptor of the outer mitochondrial membrane that is part of the translocase of the outer membrane complex. This protein is involved in the import of mitochondrial precursor proteins. [provided by RefSeq, Oct 2011]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TOMM70A translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	TOMM70A
Synonyms	TOMM70A; translocase of outer mitochondrial membrane 70 homolog A (S. cerevisiae); mitochondrial import receptor subunit TOM70; translocase of outer membrane 70 kDa subunit; mitochondrial precursor proteins import receptor
Entrez Gene ID	9868
mRNA Refseq	NM_014820

Protein Refseq	NP_055635
UniProt ID	O94826
Pathway	Metabolism of proteins; Mitochondrial protein import
Function	protein binding; protein transmembrane transporter activity