



FAIM3 blocking peptide (CDBP6346)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Fc receptors specifically bind to the Fc region of immunoglobulins (Igs) to mediate the unique functions of each Ig class. FAIM3 encodes an Fc receptor for IgM (see MIM 147020) (Kubagawa et al., 2009 [PubMed 19858324]; Shima et al., 2010 [PubMed 20042454]).[supplied by OMIM, Jul 2010]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAIM3 Fas apoptotic inhibitory molecule 3 [Homo sapiens (human)]
Official Symbol	FAIM3
Synonyms	FAIM3; Fas apoptotic inhibitory molecule 3; FCMR; TOSO; fas apoptotic inhibitory molecule 3; Fc mu receptor; IgM Fc receptor; immunoglobulin mu Fc receptor; regulator of Fas-induced apoptosis Toso
Entrez Gene ID	9214
mRNA Refseq	NM_001142473

Protein Refseq [NP_001135945](#)

UniProt ID O60667
