



TEP1 blocking peptide (CDBP6348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene product is a component of the ribonucleoprotein complex responsible for telomerase activity which catalyzes the addition of new telomeres on the chromosome ends. The telomerase-associated proteins are conserved from ciliates to humans. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TEP1 telomerase-associated protein 1 [Homo sapiens (human)]
Official Symbol	TEP1
Synonyms	TEP1; telomerase-associated protein 1; TP1; TLP1; p240; TROVE1; VAULT2; telomerase protein component 1; telomerase protein 1; p80 telomerase homolog; TROVE domain family, member 1
Entrez Gene ID	7011
mRNA Refseq	NM_007110

Protein Refseq	NP_009041
UniProt ID	Q99973
Function	ATP binding; RNA binding; contributes_to telomerase activity