



TRAP blocking peptide (CDBP6354)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TRAP triiodothyronine receptor auxiliary protein [Homo sapiens (human)]
Official Symbol	TRAP
Synonyms	TRAP; triiodothyronine receptor auxiliary protein
Entrez Gene ID	100187907