



Trim30a blocking peptide (CDBP6363)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Trim30a tripartite motif-containing 30A [Mus musculus (house mouse)]
Official Symbol	Trim30a
Synonyms	Trim30a; tripartite motif-containing 30A; Rpt1; Rpt-1; Trim30; tripartite motif-containing protein 30A; tripartite motif protein 30; tripartite motif-containing 30; regulatory protein, T lymphocyte 1; tripartite motif-containing protein 30; down regulatory protein of interleukin-2 receptor
Entrez Gene ID	20128
mRNA Refseq	NM_009099
Protein Refseq	NP_033125
UniProt ID	P15533

Function

DNA binding; metal ion binding; zinc ion binding
