



TRIM71 blocking peptide (CDBP6365)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TRIM71 tripartite motif containing 71, E3 ubiquitin protein ligase [Homo sapiens (human)]
Official Symbol	TRIM71
Synonyms	TRIM71; tripartite motif containing 71, E3 ubiquitin protein ligase; LIN41; LIN-41; E3 ubiquitin-protein ligase TRIM71; protein lin-41 homolog; abnormal cell LINEage LIN-41; homolog of C. elegans Lin-41; tripartite motif-containing 71; tripartite motif-containing protein 71
Entrez Gene ID	131405
mRNA Refseq	NM_001039111
Protein Refseq	NP_001034200
UniProt ID	Q2Q1W2

Pathway	MicroRNAs in cancer
Function	ligase activity; miRNA binding; ubiquitin-protein transferase activity; zinc ion binding