



TSHZ2 blocking peptide (CDBP6378)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TSHZ2 teashirt zinc finger homeobox 2 [Homo sapiens (human)]
Official Symbol	TSHZ2
Synonyms	TSHZ2; teashirt zinc finger homeobox 2; TSH2; ZABC2; ZNF218; OVC10-2; C20orf17; teashirt homolog 2; zinc finger protein 218; teashirt family zinc finger 2; cell growth-inhibiting protein 7; ovarian cancer-related protein 10-2; serologically defined colon cancer antigen 33 like
Entrez Gene ID	128553
mRNA Refseq	NM_001193421
Protein Refseq	NP_001180350
UniProt ID	Q9NRE2

Function	metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity
-----------------	---
