



TSHZ3 blocking peptide (CDBP6379)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TSHZ3 teashirt zinc finger homeobox 3 [Homo sapiens (human)]
Official Symbol	TSHZ3
Synonyms	TSHZ3; teashirt zinc finger homeobox 3; TSH3; ZNF537; teashirt homolog 3; teashirt 3; zinc finger protein 537; teashirt family zinc finger 3
Entrez Gene ID	57616
mRNA Refseq	NM_020856
Protein Refseq	NP_065907
UniProt ID	Q63HK5
Function	chromatin binding; metal ion binding; protein binding; sequence-specific DNA binding;

