



Tslp blocking peptide (CDBP6382)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Tslp thymic stromal lymphopoietin [Mus musculus (house mouse)]
Official Symbol	Tslp
Synonyms	Tslp; thymic stromal lymphopoietin; thymic stroma-derived lymphopoietin
Entrez Gene ID	53603
mRNA Refseq	NM_021367
Protein Refseq	NP_067342
UniProt ID	Q9JIE6
Pathway	Cytokine-cytokine receptor interaction; Jak-STAT signaling pathway

Function	cytokine activity
----------	-------------------
