



Crlf2 blocking peptide (CDBP6384)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Crlf2 cytokine receptor-like factor 2 [Mus musculus (house mouse)]
Official Symbol	Crlf2
Synonyms	Crlf2; cytokine receptor-like factor 2; CRLM2; Ly114; Tpte2; Tslpr; CRLM-2; delta 1; TSLP receptor; lymphocyte antigen 114; type I cytokine receptor delta 1; cytokine receptor-like molecule 2; thymic stromal lymphopoietin protein receptor; thymic stromal-derived lymphopoietin, receptor; transmembrane phosphoinositide 3-phosphatase and tensin homolog 2
Entrez Gene ID	57914
mRNA Refseq	NM_001164735
Protein Refseq	NP_001158207

UniProt ID	Q8CII9
Pathway	Cytokine-cytokine receptor interaction; Jak-STAT signaling pathway
Function	cytokine activity