



TTC5 blocking peptide (CDBP6391)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	TTC5 tetratricopeptide repeat domain 5 [Homo sapiens (human)]
Official Symbol	TTC5
Synonyms	TTC5; tetratricopeptide repeat domain 5; Strap; tetratricopeptide repeat protein 5; TPR repeat protein 5; stress-responsive activator of p300
Entrez Gene ID	91875
mRNA Refseq	NM_138376
Protein Refseq	NP_612385
UniProt ID	Q8N0Z6
Pathway	p53 pathway

FunctionDNA binding; chromatin binding
