



GID8 blocking peptide (CDBP6393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	GID8 GID complex subunit 8 [Homo sapiens (human)]
Official Symbol	GID8
Synonyms	GID8; GID complex subunit 8; TWA1; C20orf11; glucose-induced degradation protein 8 homolog; protein C20orf11; GID complex subunit 8 homolog; two hybrid-associated protein 1 with RanBPM
Entrez Gene ID	54994
mRNA Refseq	NM_017896
Protein Refseq	NP_060366
UniProt ID	Q9NWU2