



UHRF1BP1 blocking peptide (CDBP6400)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	UHRF1BP1 UHRF1 binding protein 1 [Homo sapiens (human)]
Official Symbol	UHRF1BP1
Synonyms	UHRF1BP1; UHRF1 binding protein 1; ICBP90; C6orf107; dJ349A12.1; UHRF1-binding protein 1; ICBP90 binding protein 1; ICBP90-binding protein 1; UHRF1 (ICBP90) binding protein 1; ubiquitin-like containing PHD and RING finger domains 1-binding protein 1
Entrez Gene ID	54887
mRNA Refseq	NM_017754
Protein Refseq	NP_060224
UniProt ID	Q6BDS2

Function

histone deacetylase binding; identical protein binding; protein binding
