



UHRF1BP1L blocking peptide (CDBP6401)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	UHRF1BP1L UHRF1 binding protein 1-like [Homo sapiens (human)]
Official Symbol	UHRF1BP1L
Synonyms	UHRF1BP1L; UHRF1 binding protein 1-like; SHIP164; UHRF1-binding protein 1-like; UHRF1 (ICBP90) binding protein 1-like; syntaxin 6Habc-interacting protein of 164 kDa
Entrez Gene ID	23074
mRNA Refseq	NM_001006947
Protein Refseq	NP_001006948
UniProt ID	A0JNW5