



ULK3 blocking peptide (CDBP6404)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ULK3 unc-51 like kinase 3 [Homo sapiens (human)]
Official Symbol	ULK3
Synonyms	ULK3; unc-51 like kinase 3; serine/threonine-protein kinase ULK3
Entrez Gene ID	25989
mRNA Refseq	NM_001099436
Protein Refseq	NP_001092906
UniProt ID	Q6PHR2
Pathway	Regulation of autophagy; TOR signaling; mTOR signaling pathway

FunctionATP binding; protein binding; protein serine/threonine kinase activity
