



ULK4 blocking peptide (CDBP6405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ULK4 unc-51 like kinase 4 [Homo sapiens (human)]
Official Symbol	ULK4
Synonyms	ULK4; unc-51 like kinase 4; FAM7C1; REC01035; serine/threonine-protein kinase ULK4; unc-51-like kinase 4
Entrez Gene ID	54986
mRNA Refseq	NM_017886
Protein Refseq	NP_060356
UniProt ID	Q96C45
Function	ATP binding; protein serine/threonine kinase activity