



UNC93A blocking peptide (CDBP6406)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	UNC93A unc-93 homolog A (C. elegans) [Homo sapiens (human)]
Official Symbol	UNC93A
Synonyms	UNC93A; unc-93 homolog A (C. elegans); Unc-93A; dJ366N23.1; dJ366N23.2; protein unc-93 homolog A; unc93 homolog A; dJ366N23.1 (putative C. elegans UNC-93 (protein 1, C46F11.1) LIKE); dJ366N23.2 (putative C. elegans UNC-93 (protein 1, C46F11.1) C-terminal LIKE)
Entrez Gene ID	54346
mRNA Refseq	NM_001143947
Protein Refseq	NP_001137419
UniProt ID	Q86WB7

Functionmolecular_function
