



Apba1 blocking peptide (CDBP6412)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Apba1 amyloid beta (A4) precursor protein binding, family A, member 1 [Mus musculus (house mouse)]
Official Symbol	Apba1
Synonyms	Apba1; amyloid beta (A4) precursor protein binding, family A, member 1; X11; Mint; Mint1; Lin-10; Urop11; X11alpha; 6430513E09Rik; amyloid beta A4 precursor protein-binding family A member 1; mint-1; UROP11-110; adapter protein X11alpha; neuron-specific X11 protein; neuronal Munc18-1-interacting protein 1
Entrez Gene ID	319924
mRNA Refseq	NM_177034
Protein Refseq	NP_796008

UniProt ID	B2RUJ5
Pathway	XPodNet - protein-protein interactions in the podocyte expanded by STRING
Function	PDZ domain binding; beta-amyloid binding; phosphatidylinositol-4,5-bisphosphate binding; protein binding; protein complex binding