



Pkdcc blocking peptide (CDBP6420)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Pkdcc protein kinase domain containing, cytoplasmic [Mus musculus (house mouse)]
Official Symbol	Pkdcc
Synonyms	Pkdcc; protein kinase domain containing, cytoplasmic; Vlk; MAD1; Adtk1; ESTM17; Sgk493; X83346; AI115348; AW548124; protein kinase domain-containing protein, cytoplasmic; sugen kinase 493; vertebrate lonesome kinase; protein kinase-like protein Sgk493
Entrez Gene ID	106522
mRNA Refseq	NM_134117
Protein Refseq	NP_598878
UniProt ID	Q5RJI4

Function

ATP binding; kinase activity; nucleotide binding; protein kinase activity; transferase activity; transferase activity, transferring phosphorus-containing groups
