



VPS53 blocking peptide (CDBP6423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein with sequence similarity to the yeast Vps53p protein. Vps53p is involved in retrograde vesicle trafficking in late Golgi. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	VPS53 vacuolar protein sorting 53 homolog (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	VPS53
Synonyms	VPS53; vacuolar protein sorting 53 homolog (S. cerevisiae); HCCS1; PCH2E; hVps53L; pp13624; vacuolar protein sorting-associated protein 53 homolog
Entrez Gene ID	55275
mRNA Refseq	NM_001128159
Protein Refseq	NP_001121631

UniProt ID	Q5VIR6
Function	protein binding