



## WIZ blocking peptide (CDBP6435)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                                                                                                                                     |
|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">WIZ widely interspaced zinc finger motifs [ Homo sapiens (human) ]</a>                                                                                                                  |
| Official Symbol | WIZ                                                                                                                                                                                                 |
| Synonyms        | WIZ; widely interspaced zinc finger motifs; ZNF803; protein Wiz; WIZ zinc finger; zinc finger protein 803; widely-interspaced zinc finger motifs; widely-interspaced zinc finger-containing protein |
| Entrez Gene ID  | <a href="#">58525</a>                                                                                                                                                                               |
| mRNA Refseq     | <a href="#">NM_021241</a>                                                                                                                                                                           |
| Protein Refseq  | <a href="#">NP_067064</a>                                                                                                                                                                           |
| UniProt ID      | O95785                                                                                                                                                                                              |

|                 |                                    |
|-----------------|------------------------------------|
| <b>Function</b> | metal ion binding; protein binding |
|-----------------|------------------------------------|

---