



YPEL2 blocking peptide (CDBP6445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	YPEL2 yippee-like 2 (Drosophila) [Homo sapiens (human)]
Official Symbol	YPEL2
Synonyms	YPEL2; yippee-like 2 (Drosophila); FKSG4; protein yippee-like 2; DiGeorge syndrome-related protein
Entrez Gene ID	388403
mRNA Refseq	NM_001005404
Protein Refseq	NP_001005404
UniProt ID	Q96QA6