



YPEL3 blocking peptide (CDBP6446)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	YPEL3 yippee-like 3 (Drosophila) [Homo sapiens (human)]
Official Symbol	YPEL3
Synonyms	YPEL3; yippee-like 3 (Drosophila); protein yippee-like 3
Entrez Gene ID	83719
mRNA Refseq	NM_001145524
Protein Refseq	NP_001138996
UniProt ID	P61236