



YPEL4 blocking peptide (CDBP6447)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	YPEL4 yippee-like 4 (Drosophila) [Homo sapiens (human)]
Official Symbol	YPEL4
Synonyms	YPEL4; yippee-like 4 (Drosophila); protein yippee-like 4
Entrez Gene ID	219539
mRNA Refseq	NM_145008
Protein Refseq	NP_659445
UniProt ID	Q96NS1